NAIPO

User Manual
3D-Rotating Shiatsu Neck
Massager with Heat

Model: MGS-150D

Web: www.naipocare.com
E-mail: support@naipocare.com
Add: Huanan City Electronic Trading Building P05 Room 116, Longgang District, Shenzhen, China

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Important Safety Instruction
When using the product, all instructions and basic precautions should always be followed, including the following:

Danger
To reduce the risk of electric shock:
Always unplug this product from the electrical outlet immediately after use and before cleaning.

Warning
To reduce the risk of burns, fire, electric shock, or injury to
01. The product should never be left unattended when plugged in. Unplug from outlet when it is not use.
02. Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
03. Do not carry this product by its supply cord or use cord as a handle.
04. Never operate the product with the air openings blocked. Keep the air openings free of lint, hair, and the like.
05. Do not drop or insert any object into any opening.
06. Do not use outdoors.
07. Do not operate where aerosol (spray) products are being used or where oxygen is administered.
08. To disconnect, switch off all controls and remove plug from outlet.
09. Close supervision is necessary when this appliance is used by, or near children, invalids, or disabled persons to ensure that they do not misuse the product.
10. Keep the cord away from heated surfaces. Do not use on sensitive skin areas or in the place of poor circulation.
11. The appliance has a heated surface. Persons who are sensitive to heat must be careful when using the appliance.

- 01 -

12. The appliance must only be supplied of safely extra low voltage corresponding to the appliance.
13. The appliance is only to be used with the power supply unit provided with the appliance.

Massage Position

- 02 -

Name and Function of The Machine Parts

1. Power button ON/OFF
2. Massage direction (Forward/Reverse)
3. Heating button: ON/OFF
4. Massage speed (Medium/High/Low)

- 03 -

How to Operate

- Connect the power adaptor to the electrical outlet.
- Press the power button to begin the massage. Heating and rolling function is fully activated.
- Press the heating button to activate heating function, press it again to stop.
- Press the direction button, the massage nodes turn anticlockwise. Press it again, the direction turns clockwise.
- Press the speed button, it will show middle / quick / slow in turn.
- Press the power button to stop the massage at any time.

To Clean and Maintain the Product
1. Cleaning
Before you clean, please make sure that the switch is in the off location. And the power plug has been pulled out from the socket. Please use dry cloth to clean the product.

2. Maintenance
Keep the produce clean and far from the wet environment.

Attention
Please do not use benzene, thinner and other aggressive detergents to clean the product or the product may fade.

Contact Us
If you need any assistance or support, please contact us. We will do our best to give you a viable solution.

Web: www.naipocare.com
Email: support@naipocare.com

Technology Parameter
Name: Neck 3D-Rotating Shiatsu Neck Massager with Heat
Model: MGS-150D
Rated Voltage: DC 12V±2.5A
Rated Power: 30 W

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